

MACEFLOR FORESTRY TRAY



8 Cells : 1500 c.c.

Tray measure: 30 x 48.5 cm.

Up cell measure: 10,2 x 10,2 cm

Bottom cell measure: 8,5 x 8,5 cm.

Height: 16 cm.

Gauge: 1,7



28 Cells: 260 c.c.

Tray measure: 30 x 50 cm.

Up cell measure: 6,5 x 6 cm.

Height: 10 cm.

Cell volume: 260 c.c.

Gauge: 1,90



28 Cells A: 350 c.c.

Tray measure: 30 x 50 cm.

Up cell measure: 6,5 x 6 cm.

Height: 13,5 cm.

Cell volume: 350 c.c.

Gauge: 1,90



35 Cells: 350 c.c.

Tray measure: 30 x 48,5 cm.

Up cell measure: 5,8 x 5,4 cm.

Height: 15,8 cm.

Cell volume: 350 c.c.

Gauge: 1,70



40 Cells: 250 c.c.

Tray measure: 30 x 50 cm.
Up cell measure: 5,5 x 5,5 cm.
Height: 11 cm.
Cell volume: 250 c.c.
Gauge: 1,90.



40 Cells A: 220 c.c.

Tray measure 30 x 48,5 cm.
Up cell measure 5,1 x 5,1 cm.
Height: 15 cm.
Cell volume: 220 c.c.
Gauge: 1,70



54 Cells: 100 c.c.

Tray measure: 30 x 50 cm.
Up cell measure: 4,8 x 4,3 cm.
Height: 8 cm.
Cell volume: 100 c.c.
Gauge: 1,90.



54 Cells A: 140 c.c.

Tray measure: 30 x 50 cm.
Up cell measure: 4,8 x 4,3 cm.
Height: 11,7 cm.
Cell volume 140 c.c.
Gauge: 1,90